exabone

Doc: EIFU02 v 6.0 (2021-06-07)

exabone

EXABONE® HA/TCP

Synthetic Bone Substitute, Granules

Synthetic Bone Substitute, Block/Disc/Wedge

EN	INSTRUCTIONS FOR USE
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
CS	NÁVOD K POUŽITÍ
DA	BRUGERVEJLEDNING
DE	GEBRAUCHSANWEISUNG
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES	INSTRUCCIONES DE USO
ET	KASUTUSJUHEND
FI	KÄYTTÖOHJEET
FR	MODE D'EMPLOI
HR	UPUTE ZA UPORABU
HU	HASZNÁLATI UTASÍTÁS
IS	NOTKUNARLEIÐBEININGAR
IT	ISTRUZIONI PER L'USO
LT	NAUDOJIMO INSTRUKCIJA
LV	LIETOŠANAS INSTRUKCIJAS
NL	GEBRUIKSAANWIJZING
NO	BRUKSANVISNING
PL	INSTRUKCJE STOSOWANIA
PT	INSTRUÇÕES DE UTILIZAÇÃO
RO	INSTRUCȚIUNI DE UTILIZARE
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
SK	NÁVOD NA POUŽITIE
SL	NAVODILA ZA UPORABO
SV	BRUKSANVISNING

CE

2797

Exabone GmbH
Route de Saint-Cergue 14, CH-1260
Nyon, Switzerland
e-mail: info@exabone.com

EC

REP

Medical Device Management Ltd
Block B, The Crescent Building, Northwood,
Santry, Dublin, D09 C6X8, Ireland

EN

INSTRUCTIONS FOR USE EXABONE® HA/TCP

Synthetic Bone Substitute, Granules
Synthetic Bone Substitute, Block/Disc/Wedge

Material

EXABONE® HA/TCP is a porous material in block form (including disc and wedge) and granular form, that is composed of nominally 60% hydroxyapatite and 40% β -tricalcium phosphate. These materials have been used successfully in clinical trials over many years. EXABONE® HA/TCP provides a synthetic alternative to human or animal cadaver bone, is free from the risk of disease transmission, and can in many cases spare the patient the trauma of autograft harvesting. EXABONE® HA/TCP is intended to slowly resorb upon implantation as part of the natural remodelling process.

Intended Purpose

EXABONE® HA/TCP synthetic bone substitute is intended to perform as an osteoconductive scaffold to support the ingrowth and fusion of adjacent viable bone when placed within an osseous defect.

Indications-for-use

EXABONE® HA/TCP is intended to be used as a scaffold to support the ingrowth of adjacent viable bone in defects that are not intrinsic to the stability of the bone structure. These defects may be located in the long bones and extremities, and EXABONE® HA/TCP granules may also be used in spine surgery (spondylolysis).

Contraindications

EXABONE® HA/TCP is contraindicated where the device is intended to provide structural support in the skeletal system and **MUST NOT BE USED TO GAIN SCREW FIXATION**. EXABONE® HA/TCP is also contraindicated in patients with:

- Acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative disease
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function
- Defects of open epiphyseal plates

Warnings

The device is provided **STERILE** by gamma irradiation, and within sealed double packaging to allow aseptic presentation. **DO NOT USE IF OPENED, PUNCTURED OR PRODUCT APPEARS DAMAGED**. Read the expiration date before use and **DO NOT USE BEYOND THE EXPIRATION DATE**. EXABONE® HA/TCP is for **SINGLE USE ONLY**. Do not attempt to re-sterilise or re-use due to the risk of cross infection.

Do not cut or reshape the implant whilst implanted in the surgical site to prevent particulate debris being introduced into the surgical site. Blocks and granules must be secured to prevent any potential migration of the implant and should only be used in procedures where the implant can be adequately contained. EXABONE® HA/TCP is opaque to x-rays and may hide areas under or above the implant on a radiograph. The device is MRI safe.

Precautions

EXABONE® HA/TCP is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement. It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support

and maintain the implant in a static, load free environment. In order to facilitate the formation of new bone, EXABONE® HA/TCP should only be implanted in direct contact with well vascularised viable bony tissue.

The device is intended for use in any patient requiring bone graft as indicated. The use of the device in paediatric patients has not been specifically evaluated however no special precautions have been identified at the time of issue. The interactions with other substances are unknown, however no special precautions have been identified at the time of issue.

Possible Complications & Adverse Reactions

Possible post-operative complications are those that may occur with any surgery. A successful result may not be achieved in every case. A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions, device displacement or device failure including fracture or crushing with or without generation of particulate debris. Allergic reaction to the product may occur in very rare circumstances.

Application

Blocks/discs/wedges may be gently shaped with a scalpel or standard surgical bone tools to adequately fit the defect and after shaping must be rinsed in sterile water or saline to remove particulate debris before attempting to implant. Prior to implantation the implant may be soaked with blood as per normal surgical practice. The device should be gently and carefully placed into the defect which should be completely filled. Secure the surgical site after implanting to prevent motion and implant displacement. If the material does not fill the site correctly, remove the implant and repeat the process with a fresh amount of EXABONE® HA/TCP

Storage

Store EXABONE® HA/TCP between 5°C and 30°C. Direct contact with sunlight or heating systems should be avoided.

Shelf Life and Disposal

The expiration date is printed on the labelling. Do not use EXABONE® HA/TCP after the expiration date. Residual material should be disposed as standard clinical waste.

Implant Card, Traceability & Event Reporting

The device is identifiable by its lot number which is pre-printed on the labelling and implant card. The implant card enclosed with this device should be supplied to the patient after implantation. Before providing it to the patient, the healthcare institution must complete the implant card with the following information:

 Name of the patient  Date of implantation

 Healthcare institution (name and address)

The summary of safety & clinical performance may be obtained via the Eudamed database <https://ec.europa.eu/tools/eudamed> where it is linked to the Basic UDI-DI: 0506017162HATCPLF, or upon request to the manufacturer.

EU: Serious incidents should be reported to the manufacturer and to the competent authority of the member state of the user and/or patient's establishment.

Other Territories: Serious incidents should be reported to the manufacturer and to the relevant national authority

Note: EXABONE® HA/TCP is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of EXABONE® HA/TCP, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.

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