 **exabone**

Doc: EIFU01 v 6.1 (2023-12-21)

# exabone

## EXABONE® PASTE

Synthetic Bone Graft, Paste

<b>EN</b>	INSTRUCTIONS FOR USE
<b>BG</b>	ИНСТРУКЦИИ ЗА УПОТРЕБА
<b>CS</b>	NÁVOD K POUŽITÍ
<b>DA</b>	BRUGERVEJLEDNING
<b>DE</b>	GEBRAUCHSANWEISUNG
<b>EL</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
<b>ES</b>	INSTRUCCIONES DE USO
<b>ET</b>	KASUTUSJUHEND
<b>FI</b>	KÄYTTÖOHJEET
<b>FR</b>	MODE D'EMPLOI
<b>HR</b>	UPUTE ZA UPORABU
<b>HU</b>	HASZNÁLATI UTASÍTÁS
<b>IS</b>	NOTKUNARLEIÐBEININGAR
<b>IT</b>	ISTRUZIONI PER L'USO
<b>LT</b>	NAUDOJIMO INSTRUKCIJA
<b>LV</b>	LIETOŠANAS INSTRUKCIJAS
<b>NL</b>	GEbruIKSAANWIJZING
<b>NO</b>	BRUKSANVISNING
<b>PL</b>	INSTRUKCJE STOSOWANIA
<b>PT</b>	INSTRUÇÕES DE UTILIZAÇÃO
<b>RO</b>	INSTRUCȚIUNI DE UTILIZARE
<b>RU</b>	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
<b>SK</b>	NÁVOD NA POUŽITIE
<b>SL</b>	NAVODILA ZA UPORABO
<b>SV</b>	BRUKSANVISNING

**CE**  
2797

Exabone GmbH  
Route de Saint-Cergue 14, CH-1260  
Nyon, Switzerland  
e-mail: [info@exabone.com](mailto:info@exabone.com)

Medical Device Management Ltd  
Block B, The Crescent Building, Northwood,  
Santry, Dublin, D09 C6X8, Ireland

**EN**

## INSTRUCTIONS FOR USE EXABONE® PASTE

Synthetic bone graft, paste

### Material

EXABONE® PASTE is a paste material composed of nominally 38% nanocrystalline hydroxyapatite in water. EXABONE® PASTE provides a synthetic alternative to human or animal cadaver bone, is free from the risk of disease transmission and can in many cases spare the patient the trauma of autograft harvesting. EXABONE® PASTE will remain as a soft paste and is intended to slowly resorb upon implantation as part of the natural remodelling process.

### Intended Purpose

EXABONE® PASTE synthetic bone graft is intended to perform as an osteoconductive scaffold to support the ingrowth and fusion of adjacent viable bone when placed within an osseous defect.

### Indications-for-use

EXABONE® PASTE is intended to be used as a scaffold to support the ingrowth of adjacent viable bone in bone defects that are not intrinsic to the stability of the bone structure. These defects may be located in the long bones and extremities, and EXABONE® PASTE may also be used for the filling of cages in spinal fusion procedures.

### Contraindications

EXABONE® PASTE is contraindicated where the device is intended to provide structural support in the skeletal system and **MUST NOT BE USED TO GAIN SCREW FIXATION**. EXABONE® PASTE is also contraindicated for use in patients with:

- Acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularised implantation site
- Uncontrolled diabetes
- Severe degenerative disease
- Collagen disease
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function
- Defects of open epiphyseal plates

### Warnings

The device is provided **STERILE** by gamma irradiation, and within sealed double packaging to allow aseptic presentation. **DO NOT USE IF OPENED, PUNCTURED OR PRODUCT OR APPLICATOR APPEARS DAMAGED OR PASTE HAS DRIED**. Read the expiration date before use and **DO NOT USE BEYOND THE EXPIRATION DATE**. EXABONE® PASTE is for **SINGLE USE ONLY**, do not attempt to re-sterilise or re-use due to the risk of cross infection. EXABONE® PASTE should only be used in procedures where the implant can be adequately contained. EXABONE® PASTE is opaque to x-rays and may hide areas under or above the implant on a radiograph. The device is MRI safe.

### Precautions

EXABONE® PASTE is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement. It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support and maintain the implant in a static, load free environment. In order to facilitate the formation of new bone, EXABONE® PASTE should only be implanted in direct contact with well vascularised bony tissue. As with any material,

care should be taken to avoid the occurrence of emboli, therefore the highly pressurised application of EXABONE® PASTE into a tightly confined defect with venous or arterial access should be avoided.

The device is intended for use in any patient requiring bone graft as indicated. The use of the device in paediatric patients has not been specifically evaluated however no special precautions have been identified at the time of issue. The interactions with other substances are unknown, however no special precautions have been identified at the time of issue.

### **Possible Complications & Adverse Reactions**

Possible post-operative complications are those that may occur with any surgery. A successful result may not be achieved in every case. A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions, device displacement or device failure. Allergic reaction to the product may occur in very rare circumstances.

### **Application**

EXABONE® PASTE is supplied as a paste in the applicator, ready for use. Optional small or large cannulas are included to assist implantation if desired. If using one of the supplied cannulas then screw the selected one to end of the syringe. The maximum dose of EXABONE® PASTE is 20cc per defect. The defect should be completely filled and the paste moulded to follow the natural contours of the bone. If using the large cannula, the remaining contents within the cannula can if necessary be applied by unscrewing and removing the syringe component and then applying the supplied pusher rod down through the cannula. Secure the surgical site after implanting to prevent motion and any implant displacement.

### **Storage**

Store EXABONE® PASTE between 5°C and 30°C, and DO NOT FREEZE. Direct contact with sunlight or heating systems should be avoided.

### **Shelf Life and Disposal**

The expiration date is printed on the label. Do not use EXABONE® PASTE after the expiration date. Residual material and the applicator should be disposed as standard clinical waste.

### **Implant Card, Traceability & Event Reporting**

The device is identifiable by its lot number which is pre-printed on the labelling and implant card. The implant card enclosed with this device should be supplied to the patient after implantation. Before providing it to the patient, the healthcare institution must complete the implant card with the following information:

 ? Name of the patient       Date of implantation  
 Healthcare institution (name and address)

The summary of safety & clinical performance may be obtained via the Eudamed database <https://ec.europa.eu/tools/eudamed> where it is linked to the Basic UDI-DI: 764015324PASTEBC, or upon request to the manufacturer.

EU: Serious incidents should be reported to the manufacturer and to the competent authority of the member state of the user and/or patient's establishment.

Other Territories: Serious incidents should be reported to the manufacturer and to the relevant national authority

*Note: EXABONE® PASTE is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of EXABONE® PASTE, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.*

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