 **exabone**

Doc: EIFU03 V8.1 (2023-12-19)

# exabone

## EXABONE® PUTTY

Synthetic Bone Graft, Putty

<b>EN</b>	INSTRUCTIONS FOR USE
<b>BG</b>	ИНСТРУКЦИИ ЗА УПОТРЕБА
<b>CS</b>	NÁVOD K POUŽITÍ
<b>DA</b>	BRUGERVEJLEDNING
<b>DE</b>	GEBRAUCHSANWEISUNG
<b>EL</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
<b>ES</b>	INSTRUCCIONES DE USO
<b>ET</b>	KASUTUSJUHEND
<b>FI</b>	KÄYTTÖOHJEET
<b>FR</b>	MODE D'EMPLOI
<b>HR</b>	UPUTE ZA UPORABU
<b>HU</b>	HASZNÁLATI UTASÍTÁS
<b>IS</b>	NOTKUNARLEIÐBEININGAR
<b>IT</b>	ISTRUZIONI PER L'USO
<b>LT</b>	NAUDOJIMO INSTRUKCIJA
<b>LV</b>	LIETOŠANAS INSTRUKCIJAS
<b>NL</b>	GEbruIKSAANWIJZING
<b>NO</b>	BRUKSANVISNING
<b>PL</b>	INSTRUKCJE STOSOWANIA
<b>PT</b>	INSTRUÇÕES DE UTILIZAÇÃO
<b>RO</b>	INSTRUCȚIUNI DE UTILIZARE
<b>RU</b>	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
<b>SK</b>	NÁVOD NA POUŽITIE
<b>SL</b>	NAVODILA ZA UPORABO
<b>SV</b>	BRUKSANVISNING

   
2797 0086

Exabone GmbH  
Route de Saint-Cergue 14, CH-1260  
Nyon, Switzerland  
e-mail: [info@exabone.com](mailto:info@exabone.com)

 EC REP

Medical Device Management Ltd  
Block B, The Crescent Building, Northwood,  
Santry, Dublin, D09 C6X8, Ireland

## INSTRUCTIONS FOR USE EXABONE® PUTTY

Synthetic bone graft, putty

### Material & Performance Characteristics

EXABONE® PUTTY is a resorbable granulated paste material composed by mass of 32% nanocrystalline hydroxyapatite 51.5% water and 16.5% HA/TCP granules (0.5 – 1.0mm). The HA/TCP granules are composed of 60% Hydroxyapatite and 40% tricalcium phosphate. The device remains as a soft paste and will perform as an osteoconductive scaffold for at least 12 months, following which the residual material will continue to slowly resorb as part of the natural remodelling process, within a period of 24 months from implantation.

### Intended Purpose

EXABONE® PUTTY is a synthetic bone graft intended to fill traumatically or surgically created bone defects that are not intrinsic to the stability of the bony structure.

### Indications-for-use

EXABONE® PUTTY is indicated for treating defects of the long bones and extremities, and EXABONE® PUTTY may also be used for the filling of cages in spinal fusion procedures.

### Clinical Benefit

EXABONE® PUTTY supports the ingrowth and fusion of adjacent viable bone when placed within an osseous defect and provides a synthetic alternative or supplement to human or animal cadaver bone that is free from the risk of disease transmission and can in many cases spare the patient the trauma of autograft harvesting.

### Contraindications

EXABONE® PUTTY is contraindicated where the device is intended to provide structural support in the skeletal system and **MUST NOT BE USED TO GAIN SCREW FIXATION**. EXABONE® PUTTY is also contraindicated for use in patients with:

- Acute or chronic infections, particularly at the implantation site
- Severe vascular or neurological disease
- Poorly vascularised implantation site
- Uncontrolled diabetes
- Severe degenerative disease
- Collagen disease
- Hypercalcemia, abnormal calcium metabolism
- Inflammatory bone disease
- Malignant tumours
- Severely impaired renal function
- Defects of open epiphyseal plates

### Warnings

The device is provided **STERILE** by gamma irradiation, and within sealed double packaging to allow aseptic presentation. **DO NOT USE IF OPENED, PUNCTURED OR PRODUCT OR APPLICATOR APPEARS DAMAGED OR PASTE HAS DRIED.** Read the expiration date before use and **DO NOT USE BEYOND THE EXPIRATION DATE**. EXABONE® PUTTY is for **SINGLE USE ONLY**, do not attempt to re-sterilise or re-use due to the risk of cross infection. EXABONE® PUTTY should only be used in procedures where the implant can be adequately contained within the defect. EXABONE® PUTTY is opaque to x-rays and may hide areas under or above the implant on a radiograph. The device is MRI safe.

### Precautions

It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support and maintain the implant in a static, load free environment. In order to facilitate the formation of new bone, EXABONE® PUTTY should only be implanted in direct contact with well vascularised bony tissue. As with any material, care should be taken to avoid the occurrence of emboli, therefore the highly pressurised application of EXABONE® PUTTY into a tightly confined defect with venous or arterial access should be avoided.

The use of the device in paediatric patients has not been specifically evaluated however no special precautions have been identified at the time of issue.

The interactions with other substances are unknown, however no special precautions have been identified at the time of issue.

## Intended User & Intended Patient

EXABONE® PUTTY is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement.

EXABONE® PUTTY is intended for use in any adult patient requiring bone graft as indicated.

## Possible Complications & Adverse Reactions

Possible post-operative complications are those that may occur with any surgery, including delayed wound healing. A successful result may not be achieved in every case. A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions, non-union, delayed or incomplete union, fracture of the newly formed bone or device displacement. Allergic reaction to the product may occur in very rare circumstances

## Application

EXABONE® PUTTY is supplied as a paste in the applicator, ready for use. Optional small or large cannulas are included to assist implantation if desired. If using one of the supplied cannulas then screw the selected one to end of the syringe. The maximum dose of EXABONE® PUTTY is 20 cubic centimetres per defect. The defect should be completely filled and the paste moulded to follow the natural contours of the bone. If using the large cannula, the remaining contents within the cannula can if necessary be applied by unscrewing and removing the syringe component and then applying the supplied pusher rod down through the cannula. Secure the surgical site after implanting to prevent motion and any implant displacement.

## Storage

Store EXABONE® PUTTY between 5°C and 30°C, and DO NOT FREEZE. Direct contact with sunlight or heating systems should be avoided.

## Shelf Life and Disposal

The expiration date is printed on the label. Do not use EXABONE® PUTTY after the expiration date. Residual material and the applicator should be disposed as standard clinical waste.

## Implant Card, Traceability & Event Reporting

The device is identifiable by its lot number which is pre-printed on the labelling and implant card. The implant card enclosed with this device should be supplied to the patient after implantation. Before providing it to the patient, the healthcare institution must complete the implant card with the following information:

Name of the patient

Date of implantation

Healthcare institution (name and address)

The summary of safety & clinical performance may be obtained via the Eudamed database <https://ec.europa.eu/tools/eudamed> where it is linked to the Basic UDI-DI: 764015324PUTTYH5, or upon request to the manufacturer.

EU: Serious incidents should be reported to the manufacturer and to the competent authority of the member state of the user and/or patient's establishment.

Other Territories: Serious incidents should be reported to the manufacturer and to the relevant national authority

*Note: EXABONE® PUTTY is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of EXABONE® PUTTY, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.*

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